

S100A13-S100 calcium binding protein A13 Gene

Reference: PTXBC000632

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein S100-A13; S100 calcium binding protein A13		
For research use only.			

GENE DESCRIPTION	S100 calcium binding protein A13
GENE ID	6284
ORIGIN SPECIES	Human
ACCESSIONS	BC000632
VECTOR	pDONR223
NCBI SYMBOL	S100A13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagcagaaccactgacagagctagaggagtcattgagaccgtggcaccacctcttcaccttgcgaaggcaggaggccgggaaggatagcctcagcgtcaacgagttcaagagctggttaccagcagttgccccatctgctcaaggatgtggcctcttctgatgagaagatg
aagagcttgatgtgaatcaggactcggagctcaagttcaatgagtactggagattgattggggagctggccaaggaaatcaggaagaagaagacctgaagatcaggaagaagtaa