

RTN3-reticulon 3 Gene

Reference: PTXBC000634

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RTN3-A1; ASYIP; HAP; NSPL2; NSPLII; reticulon-3; ASY interacting protein; NSP-like protein 2; NSP-like protein II; homolog of ASY protein; neuroendocrine-specific protein-like 2; neuroendocrine-specific protein-like II; reticulon 3		

For research use only.

GENE DESCRIPTION	reticulon 3
GENE ID	10313
ORIGIN SPECIES	Human
ACCESSIONS	BC000634
VECTOR	pDONR223
NCBI SYMBOL	RTN3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggagccgtcggcgccactcagtcctccatctcctcgtcgtcctcggagccgagccgtccgcgcggcgccggcgaggagccagcctgccccgccctggggacgaagagctgcagctcctcgtcgggtgcacgatctgattttctggagagatgtgaagaag
actgggtttgtcttggcaccacgtgatcatgctgcttccctggcagctttcagtgatcagtggtgtttcttacctcatcctggctctctctctgtcaccatcagcttcaggatctacaagtcctcatccaagctgtacagaagtcagaagaaggccatcattcaagcctacctggacgta
gacattactctgtctcagaagctttccataaattacatgaatgctgccatgggtgcacatcaacaggggccctgaaactcattatcgtctcttctggtagaagatctggttgactcctggaagctggctgtcttcatgtggctgatgacctatgttggtgctgttttaacggaatcaccccttcaattcttg
ctgaactgctcattttcagtgctccgattgtctatgagaagtacaagaccagattgatcactatgttggcatcgcccgagatcagaccaagtcattgttgaaaagatccaagcaaaactccctggaatcgccaaaaaaagcagaataa