



Reference: PTXBC000637

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CGI-86; SDR34C1; retDSR4; retSDR4; dehydrogenase/reductase SDR family member 7; dehydrogenase/reductase (SDR family) member 7; retinal short-chain dehydrogenase/reductase 4; short chain dehydrogenase/reductase family 34C member 1; dehydrogenase/reductase		

For research use only.

GENE DESCRIPTION	dehydrogenase/reductase (SDR family) member 7
GENE ID	51635
ORIGIN SPECIES	Human
ACCESSIONS	BC000637
VECTOR	pDONR223
NCBI SYMBOL	DHRS7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]