

LEPROTL1-leptin receptor overlapping transcript-like 1 Gene

Reference: PTXBC000642

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC112; Vps55; my047; leptin receptor overlapping transcript-like 1; endospanin-2		
For research use only.			

GENE DESCRIPTION	leptin receptor overlapping transcript-like 1
GENE ID	23484
ORIGIN SPECIES	Human
ACCESSIONS	BC000642
VECTOR	pDONR223
NCBI SYMBOL	LEPROTL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaggcatcaagcttggattgttcttggaggagcaatcgactgatgtttgatgttgatgtgccctccaatacaacaatactggccctcttggctattttttacatctttcacctattccatactgcatagaagaagattagtgatgatacagatgctatgagtaacg
cttgaagggaactggccatctttctacaacgggcatgtcgtgtcagctttggactccctattgtatttgcagagcacatctgattgagtgaggagctgtgcactgtttctcacaggaaacacagtcatttgcactatactaggcttttcttggcttttggagcaatgacgactcagctggc
agcagtggtga