

CD52-CD52 molecule Gene

Reference: PTXBC000644

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD52 molecule; CD52 antigen (CAMPATH-1 antigen); CDW52; CAMPATH-1 antigen; CDW52 antigen (CAMPATH-1 antigen); HEL-S-171mP; cambridge pathology 1 antigen; epididymal secretory protein E5; epididymis secretory sperm binding protein Li 171mP; he5; human epid		

For research use only.

GENE DESCRIPTION	CD52 molecule
GENE ID	1043
ORIGIN SPECIES	Human
ACCESSIONS	BC000644
VECTOR	pDONR223
NCBI SYMBOL	CD52
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagcgttcctcttctctactcaccatcagcctcctgttatggtagacagatacaaaactggactctcaggacaaaacgacaccagcccaaccagcagccctcagcatcagcagcatgagcggaggcatttctttctctgtggccaatgccataatccacctcttctgcttcagttga