

RAD51C-RAD51 homolog C (*S. cerevisiae*) Gene

Reference: PTXBC000667

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BROVCA3; FANCO; R51H3; RAD51L2; DNA repair protein RAD51 homolog 3; RAD51-like protein 2; yeast RAD51 homolog 3; RAD51 paralog C		

For research use only.

GENE DESCRIPTION	RAD51 homolog C (<i>S. cerevisiae</i>)
GENE ID	5889
ORIGIN SPECIES	Human
ACCESSIONS	BC000667
VECTOR	pDONR223
NCBI SYMBOL	RAD51C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcgcgggaagacgtccgtttgaaatgcagcgggatttggtgagtttccgctgtctccagcgggtcgggtgaagctggtgtctgcgggttcagactgctgaggaaactctagaggtgaaacccctcgagcttagcaagaagttgggatatctaaagcagaagccttagaaactct
gcaaattatcagaagaatgtctcacaataaaccaagatgctgtacatctgagtcacacaagaagtgtacagcactggaactctgtgagcaggagcataccagggttcataatcacctctgttcagcactagatgatattctgggggtgagtgcccttaataaaaaaacagaa
atttgtgtgcaccagggttggaataaacacaattatga