

### TSPAN3-tetraspanin 3 Gene

Reference: PTXBC000704

|                |                                                                                                                                             |                     |                |
|----------------|---------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                         | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                        | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | TM4-A; TM4SF8; TSPAN-3; tetraspanin-3; tetraspan 3; tetraspan TM4SF; tetraspanin TM4-A; transmembrane 4 superfamily member 8; tetraspanin 3 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | tetraspanin 3                                                                                     |
| <b>GENE ID</b>          | 10099                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC000704                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | TSPAN3                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

#### ORF (Only for reference, product provided without STOP codon )

atgggccagtgccgcatcacctctccaagaccgtgctgtcttctcaactcatcttctggggggcagctggcattttatgctatgtgggagcctatgtcttcacacttatgatgactatgaccactctttgaagatgtgtacagctcatccctgctgtagtcatagctgtaggagccctg  
cttttcatcattgggctaattggctgctgtgccaatccgggaaagtcgctgtggacttgcctttgtcatcctgctcttggttttgtcacagaagttgtgtagtgggtttgggatgtttacagagcaaaagtggaataagaggttgatcgagcattcagaaggtgtataagacctacaatg  
gaaccaacctgatgctgtagccgggctattgattgtacagagacagctgattgttgggaattcacaactactcagactgggaaaatacagattgttcaaaagaaacaaaaccagagtgctccctcttagctgctgcagagagactgccagcaattgtaatggcagcctggccacc  
cttcgacctatgctgaggggtgtgaggtctagttgtgaagaagctacaagaatcatgatgatgtgatctggggcgcactggcatttgcagctattcagctgctgggcatgctgtgtgcttgcacgtgtgtgtgcagaaggagtagagatccttctacgagctcctcatcactggcgga  
acctatgcatag