

C6orf48-chromosome 6 open reading frame 48 Gene

Reference: PTXBC000706

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	D6S57; protein G8; chromosome 6 open reading frame 48		

For research use only.

GENE DESCRIPTION	chromosome 6 open reading frame 48
GENE ID	50854
ORIGIN SPECIES	Human
ACCESSIONS	BC000706
VECTOR	pDONR223
NCBI SYMBOL	C6orf48
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagagttgtaatggcccgactgttgagtgaggggagcaggggatcccaacggcttgcgtgcctttgcgcagcagccggcgccgacgtcgccgctggctgggtag