

TIMM8B-translocase of inner mitochondrial membrane 8 homolog B (yeast) Gene

Reference: PTXBC000711

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DDP2; TIM8B; mitochondrial import inner membrane translocase subunit Tim8 B; DDP-like protein; deafness dystonia protein 2; translocase of inner mitochondrial membrane 8 homolog B		

For research use only.

GENE DESCRIPTION	translocase of inner mitochondrial membrane 8 homolog B (yeast)
GENE ID	26521
ORIGIN SPECIES	Human
ACCESSIONS	BC000711
VECTOR	pDONR223
NCBI SYMBOL	TIMM8B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagttatgttggataaatgtgtggagaagccaggggaatgcctagactctgcactgaaaattgtctctccagctgtgtgacacgcttcattgacaccactcttccatcaccagtcggtttgccagattgtacagaaaggaggcagtag