

SYNGR1-synaptogyrin 1 Gene

Reference: PTXBC000731

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms:	synaptogyrin-1; synaptogyrin 1		

For research use only.

GENE DESCRIPTION	synaptogyrin 1
GENE ID	9145
ORIGIN SPECIES	Human
ACCESSIONS	BC000731
VECTOR	pDONR223
NCBI SYMBOL	SYNGR1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaaggggtgcgtacggagcgggcaaaagccggggccttcgaccctacacctgggtccggcagccgacaccatctgcgctgtgtttggctgttccatagtgtgttcggctccatcgtgaacgaggctacctcaacagcgcctccgaggggaggagtctgca
tctacaaccgcaacccaacgctgcagctatggcgtggcgtggcgtgctccttctcactgctgctgtacctggcctggacgtgtactcccgagatcagcagcgtcaaggaccgcaagaaagccgtctccgacatgggtctcggccttctggtctggttc
gtgggattctgtacctggccaaccagtggtcaggtctccaagcccaaggacaaccaactgaacgaaggacggacgacgagccggccgacctgccttctcttccatctcactggagcctgaccgagccctggccgtcggagattcaaggacctaagcttcaggagga
gtacagcacactgttccctgcctcggcacagccgtag