

RGS4-regulator of G-protein signaling 4 Gene

Reference: PTXBC000737

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RGP4; SCZD9; regulator of G-protein signaling 4; schizophrenia disorder 9		

For research use only.

GENE DESCRIPTION	regulator of G-protein signaling 4
GENE ID	5999
ORIGIN SPECIES	Human
ACCESSIONS	BC000737
VECTOR	pDONR223
NCBI SYMBOL	RGS4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgcaaaaggctgcaggctgccggcttctgcttgaggagtgcaaaagatatgaacatcggttaggtttcctgctgcaaaatctgattcctgtgaacacaattctccacaacaagaaggacaaagtgttatttccagagagtgagccaagagggaagtcaagaaatgggctgaa
tcactggaaaacctgattagtcgaatgtggctggcagcttcaaaagcttcttgaagctggaatatagtgaggagaatatgacttctgagtcagctgtggaagagtacaagaaaatcaaatcaccatctaaactaagtcccaaggccaaaagatctataatgaattcatctcagtcaggca
accaagaggtgaacctggattcttgcaccagggaagagacaagccggaacatgctagagcctacaataacctgctttagtgaggccagaagaagatttcaacctgatggagaaggattcctaccgccgttccctcaagtctcgattctatcttatttggtaaccgtccagctgtggg
gcagaaaagcagaaaggagccaagagttcagcagactgtgttccctgggtccctcagtgctcctaa