



Reference: PTXBC000740

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CCK-B; CCK2R; GASR; gastrin/cholecystokinin type B receptor; CCK-B receptor; CCK-BR; CCK2 receptor; CCK2-R; cholecystokinin-2 receptor; gastrin receptor; cholecystokinin B receptor		

For research use only.

GENE DESCRIPTION	cholecystokinin B receptor
GENE ID	887
ORIGIN SPECIES	Human
ACCESSIONS	BC000740
VECTOR	pDONR223
NCBI SYMBOL	CCKBR
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]