

MXD3-MAX dimerization protein 3 Gene

Reference: PTXBC000745

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BHLHC13; MAD3; MYX; max dimerization protein 3; Max-associated protein 3; Max-interacting transcriptional repressor MAD3; class C basic helix-loop-helix protein 13; max dimerizer 3		

For research use only.

GENE DESCRIPTION	MAX dimerization protein 3
GENE ID	83463
ORIGIN SPECIES	Human
ACCESSIONS	BC000745
VECTOR	pDONR223
NCBI SYMBOL	MXD3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaaccttgccagcaacatccaggtcctgctcagggcgccgagttcctggagcgccgtgagagagagccgagcatggttatgctccctgtgcccgcagtcgagtcaggccccatccacaggaggaagaagcgacccccaggctcctggcgcgaggacagcgggc
ggtcagtgcaaatgaactggagaagcgagggagggccagttgaagcggctgagcggtgaagcagcagatgccccctggggccgactgtgcccgtacaccacgctgagcctgctgcgcctgcccagatgacatccagaagctggaggatcaggagcagcgggc
cgacagctcaaggagaggtgctgcgcagcaagcagcagagcctgcagcgccagctggagcagctccggggctggcagggcgccgagcgggagcggctgcccggcgacagctggactcctcaggcctctctctgagcgctcagactcagaccaagaggagctggaggtgg
atgtggagagcctggtgttgggggtgagggcgagctgctgcccggcctctgctccggccagggagcacagctactgcacggcgccggcgctggtatga