

AKIRIN2-akirin 2 Gene

Reference: PTXBC000764

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf166; FBI1; dJ486L4.2; akirin-2; fourteen-three-three beta interactant 1; akirin 2		
For research use only.			

GENE DESCRIPTION	akirin 2
GENE ID	55122
ORIGIN SPECIES	Human
ACCESSIONS	BC000764
VECTOR	pDONR223
NCBI SYMBOL	AKIRIN2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgtcggagcactctgaaaaggactctgatttcgacccgtgtgagccggcggtcccggaagcgcagcgatgtgcgccattgtcggcgccacctcggcggtgcctccccgtgtcggcgccgcccaccgcccctcttcggtgcggccgctcgcgcag
aagtatctccgaatggagccatcccccttcggcgagctctcctccgctcaccacagaacaaattctgtacaacataaaacaagagtataaacgaatgcagaagagaagacattagaacgagtttccacagacagatccgtgtgtacttctgatgcacagccacatgcatttctcctca
gtggaccagcttcaccagggacttcatctgcagcatctcaccattaaaaaagaacagcccttatttactctacggcaggttgggatgatctgtgaacgtttgttgaagaacgtgaagagaaggttcgagaagaatgaagaataatgaacacaaaactgcagaacaatgatgcgttt
gtgaagtttacgatgaataatgcgacgataatggagacagcctgctagctatgtttcatga