

SNX11-sorting nexin 11 Gene

Reference: PTXBC000768

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	sorting nexin-11; sorting nexin 11		
For research use only.			

GENE DESCRIPTION	sorting nexin 11
GENE ID	29916
ORIGIN SPECIES	Human
ACCESSIONS	BC000768
VECTOR	pDONR223
NCBI SYMBOL	SNX11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggcttttggtaggatgtcggagaaccaagaacaggaggagtgattacagtcgctgttcaggacccccgagtcagaatgagggctcctggaaactcttatgtggattataagatatctccataccaacagcaaaagcctttactgccaagactcctgtgtgctggcgctaccgt
gagttcgtgtgctgagaaagcagctacagagaatgtggtttggtgctgttctgaacttctgggaagtcaaccttctcgccacctcagatgagttcattgagaagcgacgacaaggtctgcagcacttcttgaagaggtcctgcagagtggtgttctctgacagaccagttgc
acctattcctgcaaaagccagctctcgtgctgagatagaagcctgtgtccaggggccgaagtaccatgactgtgtctgatgccattctcgatatgctatgctaaactgtggctgggcccaggagagagggcagagctcttctacctggctaaaggagaccagcctaagagttgctgttct
tccaagatcgggtaggaggagctctccctcaccgcctccagtgaaagaaaggaccatttagaagtggtggctccagttgttgactctgaggttcttcttggaaagtcacactctccaccctctctcaccattatgctgtgattttggaagacccaaagagggaacctccactctcagt
ctgtgaggagggtgtggaggagatcatgctgtgcctttggacctggtcagttagaacagttttgaaaagtga