

PRPF38A-PRP38 pre-mRNA processing factor 38 (yeast) domain containing A Gene

Reference: PTXBC000777

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PRP38A; Prp38; pre-mRNA-splicing factor 38A; PRP38 pre-mRNA processing factor 38 domain containing A; pre-mRNA processing factor 38A		

For research use only.

GENE DESCRIPTION	PRP38 pre-mRNA processing factor 38 (yeast) domain containing A
GENE ID	84950
ORIGIN SPECIES	Human
ACCESSIONS	BC000777
VECTOR	pDONR223
NCBI SYMBOL	PRPF38A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatgatgtggagtcagtgaaaggagaaggagggatgagaagtggaaagagtgccatcacctgaccgcccggagaagctaccgagacttgacaagccccgtcgtctccacactgcgtacaggaggagtaggagccggtctccagaaggcggagtcgatctcc
caaaaggagaaagccccctccctcgcgagaaaggcatcgagcaagagtgccaagacgtcaccgcagcagggtcccgagatcggcggcacagatcccgttccaagtccccaggtcatcacctagtcacagacacaggagccactcaaggtctccgaaagggtctaagaagagccac
aagaagagccggagaggggaatgagtaa