

SCAND1-SCAN domain containing 1 Gene

Reference: PTXBC000785

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAZ1; SDP1; SCAN domain-containing protein 1; SCAN-related protein RAZ1; SCAN domain containing 1		

For research use only.

GENE DESCRIPTION	SCAN domain containing 1
GENE ID	51282
ORIGIN SPECIES	Human
ACCESSIONS	BC000785
VECTOR	pDONR223
NCBI SYMBOL	SCAND1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggtacgagccgagcttggcgccactgggagtcgccggcggtgccaccggagaaactggaaggagccggttcgagctcagccctgagcgttaactgtgtgggtcctcgtgccagaggcctcaccgctgccctgagccttccagtcaccaacgccggtgcctga
agccatccctacgccccgagctggcgctccgcccctggagctgcctctcgggcccgcacccgtgagcgttagcgcctcagccgaagctgaagcgcgtccacaccaggccccgcccgtctctagactcgggtcccgagacgttccgccagcgtttccggcagttccgctaccagg
atggcggtggtccccgggagggcttccggcagctgcgggagctgtcccgccagtggtgctgcggcctgacatccgcaccaaggagcagatcgtggagatgctggtgcaagagcagctgtcgcctcctgccgagcgggtcgggccccggcggatccgcgcgcacggatgtgc
gcatcactggctga