

RBM26-RNA binding motif protein 26 Gene

Reference: PTXBC000791

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ARRS2; C13orf10; PPP1R132; PRO1777; SE70-2; ZC3H17; RNA-binding protein 26; CTCL tumor antigen se70-2; acidic rich RS domain containing 2; cutaneous T-cell lymphoma tumor antigen se70-2; protein phosphatase 1, regulatory 132; RNA binding motif protein 26		

For research use only.

GENE DESCRIPTION	RNA binding motif protein 26
GENE ID	64062
ORIGIN SPECIES	Human
ACCESSIONS	BC000791
VECTOR	pDONR223
NCBI SYMBOL	RBM26
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgaggcagggggcagggcgaggcgaggtgtgcctggatcgtgtggatcacccagggcattggagattctgcattacgagagcgatagagaagatcttctcctcatttgcgcaatatggtgaaattgaagattgtcagattgatgattcctcacttcagtcagt
aa