

POLR2K-polymerase (RNA) II (DNA directed) polypeptide K, 7.0kDa Gene

Reference: PTXBC000806

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ABC10-alpha; RPABC4; RPB10alpha; RPB12; RPB7.0; hRPB7.0; hsRPB10a; DNA-directed RNA polymerases I, II, and III subunit RPABC4; DNA directed RNA polymerases I, II, and III 7.0 kda polypeptide; DNA-directed RNA polymerase II subunit K; RNA polymerase II 7.0		

For research use only.

GENE DESCRIPTION	polymerase (RNA) II (DNA directed) polypeptide K, 7.0kDa
GENE ID	5440
ORIGIN SPECIES	Human
ACCESSIONS	BC000806
VECTOR	pDONR223
NCBI SYMBOL	POLR2K
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacaccagaaggacgttcaacctccaagcagcaaccaatgatataatctgtggagagtgacacagaaaatgaaataaaatctagggatccaatcagatgcagagaatgtggatacagaataatgtacaagaaaaggactaaaagattggtcgttttgatgctcgatga