

BRP44L-brain protein 44-like Gene

Reference: PTXBC000810

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BRP44L; CGI-129; MPYCD; dJ68L15.3; mitochondrial pyruvate carrier 1; HSPC040 protein; brain protein 44-like protein		

For research use only.

GENE DESCRIPTION	brain protein 44-like
GENE ID	51660
ORIGIN SPECIES	Human
ACCESSIONS	BC000810
VECTOR	pDONR223
NCBI SYMBOL	BRP44L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgggcgcgttggtgcgaaagcggcgactatgccgaagcaaggattccgggactacatcatgagtagcactctctggggccagtagccaactgggtctccattgtgccatcaatgatgaaaaagctccagagattacagtgggcggaatgacattgccctctgtg
ctattctttgacattcatgagattgcctacaaggtagcctcgggaactggctctgttgcagccacgcaacaatgaagtagccagctcatccaggaggcggccttatcaaacacgagatgactaaacggcatctgcataa