

NDUFA5-NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 5, 13kDa Gene

Reference: PTXBC000813

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	B13; CI-13KD-B; CI-13kB; NUFM; UQOR13; NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 5; NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 5; NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 5, 13kDa; NADH-ubiquinone oxidoreductase 13 kDa		

For research use only.

GENE DESCRIPTION	NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 5, 13kDa
GENE ID	4698
ORIGIN SPECIES	Human
ACCESSIONS	BC000813
VECTOR	pDONR223
NCBI SYMBOL	NDUFA5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgggtgtgctgaagaagaccactggcctgtgggattggctgtgtgcaatactcctcacgagagcctaagaattgtacacaagattcttgatgtcttgaggaaatccctaaaaatgcagcatatagaaagtatacagaacagattacaatgagaagctggctatggtaaagcgg
aaccagatgttaaaaaattagaagaccaactcaaggcggctcaattagaagaggtgattcttcaggctgaacatgaactaaatctggcaagaaaatgagggaatggaactatggagccattagtgaagagcctcctgccgatcagtggaatggccaatataa