

MGC5590-hypothetical protein MGC5590 Gene

Reference: PTXBC000868

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein MGC5590
GENE ID	79024
ORIGIN SPECIES	Human
ACCESSIONS	BC000868
VECTOR	pDONR223
NCBI SYMBOL	MGC5590
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaagctggggagaggattgatgccagccagctgcctcacaggggtgctggaacacgtggccatgccattagcatcctgttcggcttctggacagtttcctgtgacacctacatagtcctcgttgatcagcaagataaaaggcagccctgatgttagtgcctcctctgatgagcca
tacgccagaatccagcagagcagaaggcaatgccacgcagaagaggaccaaagtcaggtgccagaagcaggctga