

ZNF174-zinc finger protein 174 Gene

Reference: PTXBC000876

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZSCAN8; zinc finger protein 174; AW-1; zinc finger and SCAN domain-containing protein 8		

For research use only.

GENE DESCRIPTION	zinc finger protein 174
GENE ID	7727
ORIGIN SPECIES	Human
ACCESSIONS	BC000876
VECTOR	pDONR223
NCBI SYMBOL	ZNF174
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagctaaatggagataactttaagctccaacactgaagcttctccaagcaagagagacacataatagccaaactagaagagaacggggccctctctgcaaaaaaactgccagatcctgagctctgccgcagagcttcagagctttgttatcaagaggtgtctggacccc
aagaggcgtctcccagctccgacagctctgccgtcagtggttgcaacccgagctgcacaccaaggagcagattttggagcttctggtgatggagcagttcctgaccatctgccccggagatccaggtcgggtcaggcatcgttccaatgagcagcaaggagattgtgacctcg
tggaagatttcacagagcatccaagaaccaaagcagtggtggcgtttgtatgcaggggcaaaaggtgctcttggagaaaactggatctcagcttggagaacaggaactgccagacttcaaccgagactcctaggagagatctcaggagagctctccagcagagcttccag
gcaggagcttatgacggctgagccccatcattgggagaaatccccactctccaagaaccaacccccaaattggctgggacagaactcttatagaaaagacagatccaaatgcccacagatgaacttccatgcaagctatggctgagtttcattgcttaa