

VPS29-vacuolar protein sorting 29 homolog (S. cerevisiae) Gene

Reference: PTXBC000880

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	VPS29, retromer complex component; vacuolar sorting protein VPS29/PEP11; DC15; DC7; PEP11; vacuolar protein sorting-associated protein 29; PEP11 homolog; hVPS29; retromer protein; vacuolar protein sorting 29 homolog; vesicle protein sorting 29; x 007 prot		

For research use only.

GENE DESCRIPTION	vacuolar protein sorting 29 homolog (S. cerevisiae)
GENE ID	51699
ORIGIN SPECIES	Human
ACCESSIONS	BC000880
VECTOR	pDONR223
NCBI SYMBOL	VPS29
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgggcacagattggtgtgattaggagatctgcacatccacacgggtgcaacagttgccagctaaattcaaaaactctgtgccaggaaaaattcagcacattctctgcacaggaaacctttgcaccaaagagagttatgactatctcaagactctggctggtgatgttcataattgtgagaggagacttcgatgagaatctgaattatccagaacagaaagttgtgactgttgacagttcaaaattggtctgatccatggacatcaagttatccatggggagatatggccagcttagccctgttcagaggcaatttgatgtggacattctatctcgggacacacacacaaaatttgaa gcaattgagcatgaaaataaattctacattaatccagggttctgccactggggcatataatgccttggaacaacaattattccatcatttgggtgatggatccaggctctacagtggtcacctatgtgtatcagctaattggagatgatgaaagtagaacgaatcgaatacaaaaaaaccttaa