

HIST1H2BK-histone cluster 1, H2bk Gene

Reference: PTXBC000893

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H2B/S; H2BFAiii; H2BFT; H2BK; histone H2B type 1-K; H2B K; H2B histone family, member T; HIRA-interacting protein 1; histone 1, H2bk; histone cluster 1, H2bk; histone family member; histone cluster 1 H2B family member k		

For research use only.

GENE DESCRIPTION	histone cluster 1, H2bk
GENE ID	85236
ORIGIN SPECIES	Human
ACCESSIONS	BC000893
VECTOR	pDONR223
NCBI SYMBOL	HIST1H2BK
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccggaaccagcgaagtccgctcccgcccaagaagggtcgaagaaagccgtgactaaggcgcagaagaaggacggcaagaagcgaagcgcagccgcaaggagagctactccgtatacgtgtacaaggtgctgaagcaggtccacccgacacggcatctcttaag
gccatgggaatcatgaactccttcgtcaacgacatctcgaacgcatcgcggtgaggtctcccgctggcgcaatacaacaagcgtctgacatcacctccaggagatccagacggcgtgctgctgctgcccgggagttggccaagcacgcccgttccgagggcaccaag
gccgtcaccaagtacaccagcctaagtaa