

RAB10-RAB10, member RAS oncogene family Gene

Reference: PTXBC000896

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB10, member RAS oncogene family; GTP-binding protein RAB10; ras-related protein Rab-10; ras-related GTP-binding protein		

For research use only.

GENE DESCRIPTION	RAB10, member RAS oncogene family
GENE ID	10890
ORIGIN SPECIES	Human
ACCESSIONS	BC000896
VECTOR	pDONR223
NCBI SYMBOL	RAB10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgaagaagacgtacgacctgctttcaagctgctcctgacggggattcggagtgagggaagacgtgcgtcttttctgttttggatgatgccttcaatactactcttttccaccataggaatagactcaagatcaaacagttgaattacaaggaaagaagatcaagctacagatat
gggatacagcagggccaggagcgtatttcacacatcacaacctcctactacagagcgcaatgggtatcatgctatgatgacatcaccaatggtaaaagtgtttgaaacatcagcaaatggcttagaacaatagatgagcatgccaatgaagatgtggaagaatgttactaggaaacaag
tgtgatattggacgacaaaagagttgtacctaaggaaggaagacagattgcaagggagcatggtattatgtttttgagactagtgcaaaagcaaatataaacatcgaaaaggcgttcctcagcttagctgaagatatcctcgaaagaccctgtaaaagagcccaacagtgaaatgta
gatatcagcagtgaggaggcgtgacaggctggaagagcaaatgctgctga