

BET1-blocked early in transport 1 homolog (*S. cerevisiae*) Gene

Reference: PTXBC000899

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Bet1 golgi vesicular membrane trafficking protein; BET1 homolog; HBET1; Bet1p homolog; Golgi vesicular membrane trafficking protein p18; blocked early in transport 1 homolog; golgi vesicular membrane-trafficking protein p18		

For research use only.

GENE DESCRIPTION	blocked early in transport 1 homolog (<i>S. cerevisiae</i>)
GENE ID	10282
ORIGIN SPECIES	Human
ACCESSIONS	BC000899
VECTOR	pDONR223
NCBI SYMBOL	BET1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggcgtgcaggcctgggtgaaggaggtacccctggaactatgggaactatgctatgctataatgtgggtatagtcctgtgaagaagaaatgagaggctcactgaaagtctgagaagcaagtaactgctataaaatctcttccattgaaataggccatgaagttaaaaccagaat
aaattattagctgaaatggattcacaatttgattccacaactggatttctaggtataaactatgggcaaaactgaagattttatccagaggagccaaacaaagctgctgtgctatatgatgctgtttctttattgtctttttatcatttattggatttaaaactgaggatga