

RAB1A-RAB1A, member RAS oncogene family Gene

Reference: PTXBC000905

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB1A, member RAS oncogene family; GTP binding protein Rab1a; RAB1; YPT1; ras-related protein Rab-1A; RAB1, member RAS oncogene family; Rab GTPase YPT1 homolog; YPT1-related protein		

For research use only.

GENE DESCRIPTION	RAB1A, member RAS oncogene family
GENE ID	5861
ORIGIN SPECIES	Human
ACCESSIONS	BC000905
VECTOR	pDONR223
NCBI SYMBOL	RAB1A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccagcatgaatcccgaatatgattatttcaagttactctgattggcgactcaggggtggaaagtcttgcctctcttaggtttgcagatgatacatatacagaagctacatcagcacaattgggtggaattcaaaataagaactatagagtttagacgggaaacaatcaagcttcaa
atatgggacacagcagggcaggaagatttcgaacaatcacctccagttattacagaggagcccatggcatcatagttgtgtatgatgtgacagatcaggagtcctcaataatgttaaacagtggctgcaggaatagatcggttatgccagtgaataatgtcaacaattgttgtaggggaaca
aatgtgatctgaccacaagaagtagtagactacacaacagcgaaggaatttgcgtattcccttggaattccgtttttggaaccagtgctaagaatgcaacgaatgtagaacagctttcatgacgatggcagctgagattaaaaagcgaatgggtcccgagcaacagctgggtgctg
agaagtccaatgttaaaatcagagcactccagtcagcagtcaggtggaggtgctgctaa