



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

PDLIM1-PDZ and LIM domain 1 Gene

Reference: PTXBC000915

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

CLIM1; CLP-36; CLP36; HEL-S-112; hCLIM1; PDZ and LIM domain protein 1; LIM domain protein CLP-36; carboxyl terminal LIM domain protein 1; elfin; epididymis secretory protein Li 112; PDZ and LIM domain 1

For research use only.

GENE DESCRIPTION	PDZ and LIM domain 1
GENE ID	9124
ORIGIN SPECIES	Human
ACCESSIONS	BC000915
VECTOR	pDONR223
NCBI SYMBOL	PDLIM1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr