



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TMUB1-transmembrane and ubiquitin-like domain containing 1 Gene

Reference: PTXBC000936

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C7orf21; DULP; SB144; transmembrane and ubiquitin-like domain-containing protein 1; dendritic cell-derived ubiquitin-like protein; hepatocyte odd protein shuttling protein; ubiquitin-like protein DULP; ubiquitin-like protein SB144; transmembrane and ubiqu		

For research use only.

GENE DESCRIPTION	transmembrane and ubiquitin-like domain containing 1
GENE ID	83590
ORIGIN SPECIES	Human
ACCESSIONS	BC000936
VECTOR	pDONR223
NCBI SYMBOL	TMUB1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacctctgattgaaggggtgggtgatgagtgaccgtctttctggtctgtctgcttctgtgtctgtccctgctgctggctcaacgacacacctgaggcggggaccaccctgccacgccctcagggaccccaacgccatccagcccgagcgacgcatctaccga
cagcatgagagggggaggccccaggggcagagacccccagctgagacacagaggtcaagctgcacagccagagccagcacgggtgtacagcaaacccgccagccccggactccccgaggagccccctgtgtctacggctgaattctccaatgattcagagcaggtggcca
ggcgctggcccacagacacattgctctcttgaagagacagatttccgcggcggaacagcaggtgctgactcatctaacagggcagctgtataggcgacgacacccagacacctggcgagcttccactctccccaactctgttctccactgcacgtgtcacgagagctggtcc
cccaaaatccccctgccgcggggctgcagccggccccctgggctggaataaggcagctctgtctccctgtgctctctgtgtgtctgtactgcagatccagtagtccggccctcttctccctgaccgccactctgggcctggccggttaccctgtctctagtc
tctctgcttctgccctgactcccgctag

More info: Please contact us for any question: contact@proteogenix.fr