

XBP1-X-box binding protein 1 Gene

Reference: PTXBC000938

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TREB-5; TREB5; XBP-1; XBP2; X-box-binding protein 1; tax-responsive element-binding protein 5; X-box binding protein 1		

For research use only.

GENE DESCRIPTION	X-box binding protein 1
GENE ID	7494
ORIGIN SPECIES	Human
ACCESSIONS	BC000938
VECTOR	pDONR223
NCBI SYMBOL	XBP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtggtggtggcagccgcgcgaacccggcggacggaccctaaagtctgtctgtcggggcagccgcctccgcggagccccggcggcagccctgcccgtcatggtgccagccagagagggggcagcccgaggcagcagcggggggctgccccagg
cgcgcaagcgacagcgccctcagcgcactgagccccgaggagaaggcgtgaggagaaactgaaaaacagagtagcagctcagactgccagagatcgaaagaaggctcgaatgagtgagctggaacagcaagtggtatgtagaagaagagaacaaaaactttgctagaaaa
tcagctttacgagagaaaaactcatgacctgtagttagaaccaggagtagaagacagcgttggggatggatgccctgtgttctgaaagaggagcggaagccaaggggaatgaagtgagggcagtgccgggtctgctgagtcgcagcactcagactacgtgcacctctgcagcag
gtgcagggccagttgtcaccctcagaacatctcccatggattctggcggtattgactctcagattcagagctgatatcctgttggcattctggacaacttgaccagctatgtttcctcaaatgccctccccagagcctgcccagcctggaggagctcccagaggtctaccagaagga
cccagttcctaccagcctcccttctctgtcagtgaggacgtcatcagccaagctggaagccattaatgaactaa