

RBCK1-RanBP-type and C3HC4-type zinc finger containing 1 Gene

Reference: PTXBC000983

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C20orf18; HOIL-1; HOIL1; PBMEI; PGBM1; RBCK2; RNF54; UBCE7IP3; XAP3; XAP4; ZRANB4; ranBP-type and C3HC4-type zinc finger-containing protein 1; HBV-associated factor 4; RBCC protein interacting with PKC1; RING finger protein 54; heme-oxidized IRP2 ubiquitin		

For research use only.

GENE DESCRIPTION	RanBP-type and C3HC4-type zinc finger containing 1
GENE ID	10616
ORIGIN SPECIES	Human
ACCESSIONS	BC000983
VECTOR	pDONR223
NCBI SYMBOL	RBCK1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggcacagccacgacgatggacgagaagaccaagaagcgtgtgggtgagcgtggagatgctcagatgcacacccgacacatctggctcacagtgcgcctgatatgacagtggcgtctcctcaaggacatgggtttctggactatggctccaccagctctgcagcagtggtg
attggcagcggcgtggcagagaccagagaccctgcactcccatgggtgctggcagaatgggacagtgccctctatctgctgtcagcccgcaacacctccctcaacccctcaggagctgcagcgggagcggcagctcggatgctggaagatctgggttcaaggacctcag
ctgcagccgctggccctctgagccaggcccccaagccgggtccccaggaacccggacggggcagccagatgcagtgccctgagccccaccgggtggcgtggcagtgccccgggtgcaccttcatcaacaagcccacgcccctggctgtgagatgtgtgcccggc
gcgccccgagggcctaccaggtccccgctcataccagcccagcagagagagcagcgcgcctggcggcgagggagggcgcgtgctcagtagcagggagtgctgcaggccacatccgcaacagcagaggcggaggtctctgcccttcattga