

COX7C-cytochrome c oxidase subunit VIIc Gene

Reference: PTXBC001005

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	cytochrome c oxidase subunit 7C, mitochondrial; cytochrome c oxidase polypeptide VIIc; cytochrome c oxidase subunit VIIc; cytochrome-c oxidase chain VIIc; cytochrome c oxidase subunit 7C		

For research use only.

GENE DESCRIPTION	cytochrome c oxidase subunit VIIc
GENE ID	1350
ORIGIN SPECIES	Human
ACCESSIONS	BC001005
VECTOR	pDONR223
NCBI SYMBOL	COX7C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttggccagagcatccggaggttcacaacctctgtgtccgtaggaccatagaggagggccctgggaagaattgccatctcagtgaaaacaagtggtcgttactagctaagatgtgttactttggatctgcatttgctacacctctctgtagtaagaccaactgcttaaa
cataa