

PDLIM3-PDZ and LIM domain 3 Gene

Reference: PTXBC001017

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ALP; PDZ and LIM domain protein 3; alpha-actinin-2-associated LIM protein; enigma homolog; PDZ and LIM domain 3		

For research use only.

GENE DESCRIPTION	PDZ and LIM domain 3
GENE ID	27295
ORIGIN SPECIES	Human
ACCESSIONS	BC001017
VECTOR	pDONR223
NCBI SYMBOL	PDLIM3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccccagacggtgatccctccggccctgcgccctgggcttcaggctctcaggggcatagacttcaaccagccttgggtcaccagcaggtacaccaggaagcaaggcggcagctgccaaacctgtgtcctggagatgtcatcctggctattgacgcttgggacagagtcctg
actcatgctgatgcgcaggacaggaattaaagcagcagctcaccagctgtgtctcaaaattgacaggggagaaactcactatgggtcccaagtagctgaagatgggaagcccatcttcaaaatcaactagaatcagaaccacaggaattcaaacccattgggtaccgcgcacaacag
aaggcccgagcctttgttcagctgcaaacattgatgacaaaagacaggtagtgagcgtctctataactcgccaattgggctctattcaactagcaatatacaagatgcgcttcacggacagctgcggggtctcattctagctcacctcaaaagacgggaactacttgaacacagcata
atattcgcccaaaccttctgga