

RPL39-ribosomal protein L39 Gene

Reference: PTXBC001019

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L39; RPL39P42; RPL39_23_1806; 60S ribosomal protein L39; ribosomal protein L39		

For research use only.

GENE DESCRIPTION	ribosomal protein L39
GENE ID	6170
ORIGIN SPECIES	Human
ACCESSIONS	BC001019
VECTOR	pDONR223
NCBI SYMBOL	RPL39
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcttcacaaagacttcaggattaagcgattcctggccaagaacaaagcaaaatcgctccattccccagtggaatcgatgaaaactggaaataaaatcaggtacaactccaaaggagacattggagaagaaccaagctgggtctataa