



Reference: PTXBC001023

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HLA-DR1B; HLA-DR3B; major histocompatibility complex, class II, DR beta 3; HLA class II histocompatibility antigen, DR beta 3 chain; MHC class II HLA-DR beta 3 chain; MHC class II antigen DR beta 3 chain; MHC class II antigen DRB3; human leucocyte antigen		

For research use only.

GENE DESCRIPTION	major histocompatibility complex, class II, DR beta 3
GENE ID	3125
ORIGIN SPECIES	Human
ACCESSIONS	BC001023
VECTOR	pDONR223
NCBI SYMBOL	HLA-DRB3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]