

EIF3K-eukaryotic translation initiation factor 3, subunit K Gene

Reference: PTXBC001031

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ARG134; EIF3-p28; EIF3S12; HSPC029; MSTP001; PLAC-24; PLAC24; PRO1474; PTD001; eukaryotic translation initiation factor 3 subunit K; eIF-3 p28; eukaryotic translation initiation factor 3, subunit 12; muscle specific; muscle-specific gene M9 protein		

For research use only.

GENE DESCRIPTION	eukaryotic translation initiation factor 3, subunit K
GENE ID	27335
ORIGIN SPECIES	Human
ACCESSIONS	BC001031
VECTOR	pDONR223
NCBI SYMBOL	EIF3K
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgatgtttgagcagatgagagccaacgtgggcaagtgtctcaagggtatcgacaggtacaatcctgagaacctggccaccctggagcgtatgtagagacgcaggccaaggaaaatgcctatgatctggaagccaacctggctgtcctgaagctgtaccagttcaaccagccttc
tttcagaccacgggtcaccgccagatcctgtgaaggccctcaccaactgcccacacagacttcacccgtgtcaagtgcatgatcgaccaggcacatcaagaagaacggccaatccgacagattttgacctcggggacctgtggagacctgccatttcaggcctctggcaagccc
tggatgaaaacatggacctcttgaaggtataactgctttgaagactctgtccgaaagtattatctgcatgtgtgggtatcacttaccagcacattgaccgctgctgctgctggccgagatgctcggggatctgtcggacagccagctaaaggtgtggaatgagcaaatacggctggagtgccg
acgagtcggggcagatcttcatctgtagccaagaagagagacattaaaccaagaacattgtggagaagattgacttgacagtgtgtccagcatcatggcctcctccagtaa