

TMEM14B-transmembrane protein 14B Gene

Reference: PTXBC001033

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 14B		
For research use only.			

GENE DESCRIPTION	transmembrane protein 14B
GENE ID	81853
ORIGIN SPECIES	Human
ACCESSIONS	BC001033
VECTOR	pDONR223
NCBI SYMBOL	TMEM14B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagaagccctctccattagtcctttgcattggtttggcttggctacacagcactgggtgtttctgtgtgggacgttggctatgtaaaacaggcagcgtgccgtccctggctgcatggctgctcttcggcagcttagccggcctgggtgcttaccagctgtatcaggatccaaggaaac
gtttggggtttcttagccgctacatctgttactttgttgggtgttatgggaatgagatcctactactatggaaaattcatgcctgtagggttaattgcagggtgccagtttgcgtatggccgcaaatgtggagttcgtatgttgatgacatctgattag