

MRPL16-mitochondrial ribosomal protein L16 Gene

Reference: PTXBC001040

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L16mt; MRP-L16; PNAS-111; 39S ribosomal protein L16, mitochondrial; mitochondrial ribosomal protein L16		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein L16
GENE ID	54948
ORIGIN SPECIES	Human
ACCESSIONS	BC001040
VECTOR	pDONR223
NCBI SYMBOL	MRPL16
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtggaggctgctgctcgctagtcgcccgtcctgctgggtgccctgtcagattcctgggcaactcctcccgccagtcgtgctgaaagacactgctccagtagtaccagtttgaagatgtttcattcctgaaaaacccaagcttagatttgaaggccaccactgtgccaaaagt
aagaagagaacctaaaaatttaagtacataggggacgtccactgaagctacggagtttacagaaggcaatttgcattctggcattgggtggtgctacactgacattggggccactttgaaatgatgcgcctgacaatcaaccgctctatggacccaagaacatgttgccatagggga
gtaccagcccttcaagcccatcactgcgaaaagtgttgggcatcgatgggggagggcaaggtgctattgaccactacgtgacacgtgtgaaggctggccgctgtgttagagatgggtggcggtgtgaattgaagaagtgaaggtttcctgaccaggttgcccaagttgcc
cttcgcagcaaaaggctgtgagccgctgggacttagaagatgcgaaaagatcaagaggaaagagaacgtacaaccagaacccctggacattgagcgaatagccactgccaacatgctgggcatacggaaagtactgagcccatgacttgaccacaaggggaatactgggg
caagttctacatgcccacacgtgttag