

TXNL4A-thioredoxin-like 4A Gene

Reference: PTXBC001046

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BMKS; DIB1; DIM1; SNRNP15; TXNL4; U5-15kD; thioredoxin-like protein 4A; DIM1 protein homolog; spliceosomal U5 snRNP-specific 15 kDa protein; thioredoxin-like 4; thioredoxin-like U5 snRNP protein U5-15kD; thioredoxin like 4A		

For research use only.

GENE DESCRIPTION	thioredoxin-like 4A
GENE ID	10907
ORIGIN SPECIES	Human
ACCESSIONS	BC001046
VECTOR	pDONR223
NCBI SYMBOL	TXNL4A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgtacatgtctccgcacctgcacaacggctggcaggtggaccaggccatctctcggaggaggaccgctgtgtcgtcatccgcttcggccacgactgggatcctacgtgcatgaagatggacgaggtcctgtacagcatgccgagaaggttaaaatttgcagttattatctgt
ggatattacagaagtgcctgacttcaacaaaatgatagtatagatccatgtactgtcatgtttttctcaggacaacacacatcatgattgacttggggactggcaacaacaagaattaactggccatggaggacaagcaggagatggtggacatcatcgagacggtgtaccgcgg
ggcccgaaaggccgctggtgtgtcctcccaaggactactccaccaagtaccgctactga