

TMEM9-transmembrane protein 9 Gene

Reference: PTXBC001106

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DERM4; TMEM9A; transmembrane protein 9; dermal papilla-derived protein 4		

For research use only.

GENE DESCRIPTION	transmembrane protein 9
GENE ID	252839
ORIGIN SPECIES	Human
ACCESSIONS	BC001106
VECTOR	pDONR223
NCBI SYMBOL	TMEM9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagctcttcttgggtgctgtggtcgggtgttctgtggtgccccagctgaagccaagaagttctgaagatatccggtgcaaatgcatctgtccacctatagaacacatcagtggtgacatttacaaccagaatgtatccagaaggactgcaactgacctgacgtgtggagccca
tgccagtgccctggccatgacgtggagggcctactgctgtgctgagtgaggtacgaggagcgcagcaccaccacatcaaggatcatctacgtgctggtgggtgacctgtgtgctctacatggccttctgatgtgtggacctctgatccgaaagccggatgcatac
actgagcaactgcacaatgaggaggagaatgaggatgctgctctatggcagcagctgtgcatccctcgggggacccgagcaaacacagtcctggagcgtgtggaagggtgccagcagcgggtggaagctgcaggtgcaggagcagcggagacagcttctgatcggcacaaga
tgctcagctag