

AP1S2-adaptor-related protein complex 1, sigma 2 subunit Gene

Reference: PTXBC001117

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DC22; MRX59; MRXS21; MRXS5; MRXS5F; PGS; SIGMA1B; AP-1 complex subunit sigma-2; adapter-related protein complex 1 sigma-1B subunit; adaptor protein complex AP-1 sigma-1B subunit; adaptor-related protein complex 1 subunit sigma-1B; clathrin adaptor complex		

For research use only.

GENE DESCRIPTION	adaptor-related protein complex 1, sigma 2 subunit
GENE ID	8905
ORIGIN SPECIES	Human
ACCESSIONS	BC001117
VECTOR	pDONR223
NCBI SYMBOL	AP1S2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagtttatgttcttttagtcgtcagggaaagcttcgactgcaaaaatggtatgtccactatcagacaaagagaagaaaaagatcacagaagaactgttcagaccgttttagcacggaaacctaaaatgtgcagctccttgagtgccgagatctgaagattgtttacaaaagatatgcta
gtctgtattttgtctgtctattgaggatcaggacaatgaactaattaccctggaataattcatcgttatgtggaattacttgacaagattttcggcagtgctgtggaactagatcatcttaatttgagaagccttattttatgttgatgatttctttggaggagggaagttcaggaaacatccaag
aaaatgtcttaagcaattgagcagcgtgatctactgcaggagggaagctgaaacccacgtatgttcttgaagaaattggactgacataa