

LGALS3-lectin, galactoside-binding, soluble, 3 Gene

Reference: PTXBC001120

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CBP35; GAL3; GALBP; GALIG; L31; LGALS2; MAC2; galectin-3; 35 kDa lectin; IgE-binding protein; MAC-2 antigen; advanced glycation end-product receptor 3; carbohydrate-binding protein 35; galactose-specific lectin 3; laminin-binding protein; lectin L-29; lec		

For research use only.

GENE DESCRIPTION	lectin, galactoside-binding, soluble, 3
GENE ID	3958
ORIGIN SPECIES	Human
ACCESSIONS	BC001120
VECTOR	pDONR223
NCBI SYMBOL	LGALS3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagacaattttcgtccatgatcgttatctgggtctggaacccaaccccaagatggcctggcgcgatggggaaccagcctgctggggcagggggctaccaggggcttctatcctggggcctacccggggcaggcaccacaggggcttctggacaggcacctcc
aggcgctaccctggagcacctggagcttatcccgagcacctgacacgtggagtctaccaggccacccagcggccctggggcctaccatctctggacagccaagtggccacggagcctacccctggccctatggcgccctgctgggcccactgattgtcctataacc
tgcctttgcctggggagtggtgctcgcacgtgataacaattctgggcacgtggaagcccaatgcaaacagaattgctttagatttccaaagagggaatgatgttgccttccactttaaccacgcttcaatgagaacaacaggagagtcattgtttgcaatacaaaagctgataataactgg
ggaagggaagaagacagtcggtttccattgaaagtgggaaccattcaatacaagtaactggtgaacctgaccactcaaggtgcagtgatgatgtcacttgttcagtaaatcatcgggttaaaaaactcaatgaatcagcaaacctgggaattctgtgacatagacctcac
cagtgcttcatataccatgatataa