

DDX54-DEAD (Asp-Glu-Ala-Asp) box polypeptide 54 Gene

Reference: PTXBC001132

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATP-dependent RNA helicase DDX54; ATP-dependent RNA helicase DP97; DEAD (Asp-Glu-Ala-Asp) box polypeptide 54; DEAD box RNA helicase 97 kDa; DEAD box helicase 97 kDa; DEAD box protein 54; DEAD-box helicase 54		

For research use only.

GENE DESCRIPTION	DEAD (Asp-Glu-Ala-Asp) box polypeptide 54
GENE ID	79039
ORIGIN SPECIES	Human
ACCESSIONS	BC001132
VECTOR	pDONR223
NCBI SYMBOL	DDX54
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcacctcggatgtggagccgacacccgggagatgggtgcgtgccagaacaagaagaagaaggtctggaggtccagtcctatgggcctgagctacccgggtgttcaaggcatcatgaagaagggggagcctttgagcagcaggcagctggcgctgtcctggacttgatggg
ggatgaagcccagaacctgacgaggggcccgcagcagctcaagtgggaccgttaagaagaagcgggtttgtgggacagtcaggacaggaagacaagaagattaaagacagagagcggccgctacatcagcagctctacaagcgagacctctatcagaagtggaaacagaaaca
gaaaattgatgctgactcggacgaagaaggggcatctgaccggcgagggccagagcgaagaggtgggaagcagacgtggccaagcaggtgcatccggcccccacgccccagccagccctgcaggccgagtcggccgggaactcaagaccaagcagcagatcctgaag
cagcggcgccgggcccagaagctgcacttctgcagcgtgggtggcctcaagcagctctctgccgcgaaccgcccgcgtccaggagctgcagcagggcgcccttcggccgggtgcccgtccaagaagggcaagatgcggaagaggtga