

DNAJB11-DnaJ (Hsp40) homolog, subfamily B, member 11 Gene

Reference: PTXBC001144

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ABBP-2; ABBP2; DJ9; Dj-9; EDJ; ERdj3; ERj3; ERj3p; PRO1080; UNQ537; dnaJ homolog subfamily B member 11; APOBEC1-binding protein 2; DnaJ (Hsp40) homolog, subfamily B, member 11; DnaJ protein 9; ER-associated DNAJ protein 3; ER-associated Hsp40 co-chaperone		

For research use only.

GENE DESCRIPTION	DnaJ (Hsp40) homolog, subfamily B, member 11
GENE ID	51726
ORIGIN SPECIES	Human
ACCESSIONS	BC001144
VECTOR	pDONR223
NCBI SYMBOL	DNAJB11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctccgcagaaactgagcaccctttgctgtgtctatccctcagggggcgggtgattgcccggagagattctataagatcttgggggtgctcgaagtgcctctataagggatataaaaggcctataggaactagccctgcagcttcacccgaccgggaacctgatgatcac
aagcccaggagaaatccaggatctgggtgctgtctatgaggttctgtcagatagtgagaaacggaaacagtagcagatactatggfagaagaaggattaaagatgggtcatcagagctcccatggagacattttcacacttctttggggtattgtttcatgtttggaggaaacctctgcagc
aagacagaaatattccaaggagagtgatattatgtagatcagaagtcactttggaagaagtatatgcaggaaattttggaagtagtagaacaacacctgtgcaagcgaggctcctgcaaacggaagtgcattgtcggcaagagatgcggaccaccagctggccctggg
cttccaaatgacccaggaggtggtctgcgacgaatgcctaatgtcaactagtagaatgaagaacgaacgctgggaagtagaataagagcctgggtgagagacggcatggagtaccccttattggagaagtgagcctcacgtggatggggagcctggagattacgttcgaatcaa
agttgtcaagcaccatatttgaaggagaggagatgattgtacacaatgtgacaatcattagttgagtcactggttggctttgagatgagatattactcacttgatggtcacaaagggtacatatcccggaataagatcaccagccaggagcgaagctatggaagaagggaagg
ctcccaactttgacaacaacaatatcaagggtctttgataatcattttgatgtgatttccaaaagaacagttacagaggaagcgagagaagggtatcaaacagctactgaacaagggtcagtcagaaggtatacaatggactgcaaggatattga