

RAB31-RAB31, member RAS oncogene family Gene

Reference: PTXBC001148

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB31, member RAS oncogene family; Rab22B; ras-related protein Rab-31; ras-related protein Rab-22B		

For research use only.

GENE DESCRIPTION	RAB31, member RAS oncogene family
GENE ID	11031
ORIGIN SPECIES	Human
ACCESSIONS	BC001148
VECTOR	pDONR223
NCBI SYMBOL	RAB31
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatggcgatacgggagctcaagtgtgcttctcgggacactgggggtgggaatcaagcatcgtgtgtcgtatgtccaggatcattgaccacaacatcagccctactattggggcatcttttatgacaaaactgtgcttggaaatgaacttcacaagtctctatctgggacact
gctggtcaggaaacgggtttcattcattggtcctcatgtactatcagggtcagctgagctgtatcgtgtatgatattaccaagcaggatcattttatacctgaagaaatgggtcaaggagctgaaagaacatggtccagaaaacattgtaatggccatcgtctgaaacaagtgcgacctctca
gatattaggaggttccccgaaggatgctaaggatcgtgaatccataggtgccatcgtgttgagacaagtgcacaaaatgctattaataatcgaagagctcttcaaggatcagccgagatccacccttgacccccatgaaaatggaacaatggaacaatcaagtgtgagaa
gccaacctgcaagccagccgccgtgtgtga