



Reference: PTXBC001149

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UTP14A small subunit processome component; UTP14A small subunit (SSU) processome component; NYCO16; SDCCAG16; dJ537K23.3; U3 small nucleolar RNA-associated protein 14 homolog A; UTP14, U3 small nucleolar ribonucleoprotein, homolog A; antigen NY-CO-16; ser		

For research use only.

GENE DESCRIPTION	UTP14, U3 small nucleolar ribonucleoprotein, homolog A (yeast)
GENE ID	10813
ORIGIN SPECIES	Human
ACCESSIONS	BC001149
VECTOR	pDONR223
NCBI SYMBOL	UTP14A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]