

RAB33A-RAB33A, member RAS oncogene family Gene

Reference: PTXBC001157

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB33A, member RAS oncogene family; RabS10; ras-related protein Rab-33A; Small GTP-binding protein S10		

For research use only.

GENE DESCRIPTION	RAB33A, member RAS oncogene family
GENE ID	9363
ORIGIN SPECIES	Human
ACCESSIONS	BC001157
VECTOR	pDONR223
NCBI SYMBOL	RAB33A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgagcccatctgggagcctgcagccgcctggcgctggcgtgcctggagctcgactcgtcgtggaccagtagtcagattcgcatctcaaaataatcgtgattgggactccaacgtgggcaagacctgcctgacctccgtctgcgggggtacc
ttccagacaagactgaagccaccatcgcgctggactcagggaagaccgtggaaatcgaggcgagaaatcaagggtcaggtgtggacacagcagggtcaggaaacgtttccgcaaaagcatggctcagcatctactaccgcaacgtacatgccgtggtcttcgtatgacgtcac
caagatgacatctttcaccaacctcaaaatgtggatccaagaatgcaatgggcatgctgtgccccactagtcaccaaggtcgttgggcaacaagtgtgacttgagggaacagatccagggtgcctccaactagccctgaaattgtgatgccacaacatgctctgtttgagacatcg
gccaaggaccccaagagagccagaacgtggagtcgattttcatgtgttggcttgcgattgaagcccagaaaatccctgctgtatcgtgatcgtgagaggcagcaggggaaggtgcagaaactggagttcccaggaagctaacagtaaaacttccttcctgttga