

U2AF1-U2 small nuclear RNA auxiliary factor 1 Gene

Reference: PTXBC001177

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FP793; RNU2AF1; U2AF35; U2AFBP; splicing factor U2AF 35 kDa subunit; U2 small nuclear RNA auxiliary factor 1; U2 small nuclear ribonucleoprotein auxiliary factor, 35-KD subunit; U2 snRNP auxiliary factor small subunit; U2(RNU2) small nuclear RNA auxiliary		

For research use only.

GENE DESCRIPTION	U2 small nuclear RNA auxiliary factor 1
GENE ID	7307
ORIGIN SPECIES	Human
ACCESSIONS	BC001177
VECTOR	pDONR223
NCBI SYMBOL	U2AF1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggagtatctggcctccatctcggcaccgagaagacaagtaactgttcattttatcctaaattggagcatgtcgtcatggagacaggtgctcgtggtgcacaataaacgacgtttagccagaccattgccctcttgaacattaccgtaaccctcaaaactctccagctctgtg
acgggttgcgctgtgccgtgagcgtatgtggagatgcaggaaactatgatgagtttttggagggtttttacagaaatggaggagaagatggggaagtagaggagatgaacgtctgtgacaacctgggagaccacctgggtgggaacctgtacgtcaagtttcgctgagggaagatgc
ggaaaaggctgtgattgacttgaataaccgttggtttaattggacagccgatccagccgagctgtcaccgtgacggacttcagagaagcctgctgctcagtagatgagggagaatgcacagaggcggcttctgcaacttcattgaagccattccagagagctcgccggg
gagctgtatggccgctgcgaagaagcatagatcaagatccgatccgggagcgtgttctcgtctagagaccgtgtgtggtggcgtggcgtggtggtgaggtggcggcgacgggagcgtgacaggagcggtcgagagatctgaaagatctggcgattctga