

MRPS34-mitochondrial ribosomal protein S34 Gene

Reference: PTXBC001182

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MRP-S12; MRP-S34; MRPS12; 28S ribosomal protein S34, mitochondrial; S34mt; mitochondrial 28S ribosomal protein S34; mitochondrial ribosomal protein S34		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein S34
GENE ID	65993
ORIGIN SPECIES	Human
ACCESSIONS	BC001182
VECTOR	pDONR223
NCBI SYMBOL	MRPS34
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgcggaagaaggtgctccgctgatcgaggctggccgcgcgtgcgcgcctcgggagcaactgaacaggccgcgcgactccagctctacgcggtggactacgagaccttgacgcggccttcttgacgccggctgccgtccggcctgggacgtg
cgccgcgagagccgcctcttgagctgctcggccctccgctcttcggcctggccgcctggtcagcgcgaagtcctgctgtggcagcagcagagccgtgctactggcctcagcgggtgcggcccgactacagggcgagaaactggaccacgggaaggcctgggca
tcctgacctcaaaagggaagactgagagcagggcgggagatcgaacacgtcatgtacatgactggcggctggtgcccgaagcagcagaggagccttcacgcggttcacgccggcgggaagacagcctggcctccgtgccgtaccgcctctcctccggccatgattatcg
cagaacgacagaaaaatggagacacacagcagcagggagccatgctgaatgtgcagaggatagcagtggaacctgggattacctgcacaaacaggaagacaaggaaggccaaggcaccctctag